



## G2P Plant DNA Purification Kit (Spin Column Based)

Cat#  PUR16-50  PUR16-250  PUR16-50-SPack Size:  50 reactions  250 reactions  5 reactions

Storage: Room Temperature

### Introduction

The G2P Plant DNA purification kit is a state-of-the-art solution that addresses the challenges of DNA extraction from plant tissue, allowing for the recovery of high-quality DNA from freshly growing plant leaves. With advanced features, this kit offers improved removal of inhibitors and increased quantities of plant and plant pathogen DNA, resulting in high-performance outcomes in downstream applications. The silica-based DNA isolation process is both rapid and straightforward, utilizing a spin column configuration. Depending on the type of sample utilized, typical yields range from 3 to 15µg of high-quality DNA.

### Kit Content

- Buffer PDC - 30mL
- RNase A (10mg/mL) - 500µL (Store @ 4°C)
- Buffer D - 10mL
- Buffer E - 5mL
- 5  50 Spin columns
- 5  50 Collection tubes

### Reagents/material required but Not Provided


- Chloroform:Isoamyl alcohol (24:1)
- 100% Ethanol
- Plastic Pestle/Tissue Lyser: For Mechanical lysis of the plant tissue

### Protocol


**Note: Add 40 ml (1:4) of absolute ethanol in Buffer D (Wash buffer) before starting the protocol**

1. Take ~50-80mg fresh growing leaf/pulp/soft stem tissue into a 1.5ml microcentrifuge tube
2. Add 600 µl of Buffer PDC and homogenize well using plastic pestle  
*Note: A good homogenization of the plant tissue is must to get a good quality and quantity of DNA. We highly recommend to use plastic pestle coupled with mechanical homogenizer/tissue lyser. If homogenizer/lyser is not available, use the liquid nitrogen base tissue pulverization method.*
3. Incubate at 55°C for 30 min  
*Note: If the tissue is hard homogenize again after 10 min of incubation*
4. Add 10 µl of RNase A and incubate at 37°C for 30 min
5. Add 600 µl of Chloroform:Isoamyl alcohol (24:1) in the above mix and invert mix
6. Centrifuge at 12000 x g for 10 min at 4°C
7. Transfer the upper aqueous layer into a fresh 1.5ml microcentrifuge tube
8. Repeat step 5-7 once again (Optional step for better quality)
9. Add 1 volume 100% chilled isopropanol invert mix and load it on the column and wait for 2 min
10. Centrifuge at 12000 x g for 5 min at 4°C  
*Note: If the entire volume has not passed through the column then increase the centrifugation time*
11. Discard the flow through which is collected in the collection tube
12. Add 500 µl of buffer D on the spin column and centrifuge at 12000 x g for 2 min at 4°C
13. Discard the flow through which is collected in the collection tube
14. Repeat step 12-13 once again
15. Dry spin the spin column by centrifuge at 12000 x g for 5 min at 4°C
16. Put the spin column into a fresh 1.5 ml microcentrifuge tube and add 50 µl of buffer E
17. Centrifuge at 10000 x g for 2 min to collect the pure DNA

## Graphical Protocol



## Result



**Fig:** Various plant genomic DNA was isolated using PUR17 (lane 1-4) and PUR16 (lane 6-9). 100ng of DNA sample from each plant was loaded on 0.8% agarose gel.