

## G2P Plant DNA Purification Kit (Spin Column Based)



Cat#  PUR16-50  PUR16-250  PUR16-50-S

Pack Size:  50 Preps,  250 Preps,  5 Preps

Storage: Room Temperature\*

### Introduction:

The G2P Plant DNA purification kit is a state-of-the-art solution that addresses the challenges of DNA extraction from plant tissue, allowing for the recovery of high-quality DNA from freshly growing plant leaves. With advanced features, this kit offers improved removal of inhibitors and increased quantities of plant and plant pathogen DNA, resulting in high-performance outcomes in downstream applications. The silica-based DNA isolation process is both rapid and straightforward, utilizing a spin column configuration. Depending on the type of sample utilized, typical yields range from 3 to 15µg of high-quality DNA.

**Kit contents:** \* Store at 4°C

S.No	Components	PUR16-50	PUR16-250	PUR16-50-S
1	Buffer PDC	30mL	150mL	3mL
2	Rnase A(10mg/mL)*	500 µl	2.5mL	50 µl
4	Buffer D (Wash Buffer)	10mL	50mL	2mL
5	Buffer E (Elution Buffer)	5mL	25mL	500 µl
6	Spin Columns and Collection Tubes	50	250	5

### Reagents/material required but Not Provided

- Chloroform:Isoamyl alcohol (24:1)
- 100% Ethanol
- Plastic Pestle/Tissue Lyser: For Mechanical lysis of the plant tissue

### Protocol

**Note: Add 40 ml (1:4) of absolute ethanol in Buffer D (Wash buffer) before starting the protocol**

1. Take ~50-80mg fresh growing leaf/pulp/soft stem tissue into a 1.5ml microcentrifuge tube
2. Add 600 µl of Buffer PDC and homogenize well using plastic pestle


**Note:** A good homogenization of the plant tissue is must to get a good quality and quantity of DNA. We highly recommend to use plastic pestle coupled with mechanical homogenizer/tissue lyser. If homogenizer/lyser is not available, use the liquid nitrogen base tissue pulverization method.

3. Incubate at 55°C for 30 min  
*Note: If the tissue is hard homogenize again after 10 min of incubation*
4. Add 10 µl of RNase A and incubate at 37°C for 30 min
5. Add 600 µl of Chloroform:Isoamyl alcohol (24:1) in the above mix and invert mix
6. Centrifuge at 12000 x g for 10 min at 4°C
7. Transfer the upper aqueous layer into a fresh 1.5ml microcentrifuge tube
8. Repeat step 5-7 once again (Optional step for better quality)
9. Add 1 volume 100% chilled isopropanol invert mix and load it on the column and wait for 2 min
10. Centrifuge at 12000 x g for 5 min at 4°C


**Note:** If the entire volume has not passed through the column then increase the centrifugation time

11. Discard the flow through which is collected in the collection tube
12. Add 500 µl of buffer D on the spin column and centrifuge at 12000 x g for 2 min at 4°C
13. Discard the flow through which is collected in the collection tube
14. Repeat step 12-13 once again
15. Dry spin the spin column by centrifuge at 12000 x g for 5 min at 4°C
16. Put the spin column into a fresh 1.5 ml microcentrifuge tube and add 50 µl of buffer E
17. Centrifuge at 10000 x g for 2 min to collect the pure DNA

## Graphical Protocol



## Result



**Fig:** Various plant genomic DNA was isolated using PUR17 (lane 1-4) and PUR16 (lane 6-9). 100ng of DNA sample from each plant was loaded on 0.8% agarose gel.