



Gene to Protein Pvt. Ltd.

Gene to Protein Pvt. Ltd.

www.genetoprotein.com

info@genetoprotein.com

800 GENOME, 800 GENETIC

Lambda DNA/HindIII Ladder

Cat# L01

Pack Size: 100 µg

Concentration: 500 mkg/ml

Storage buffer (TE): 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA.

Storage: -20°C.

Usage recommendation: Dilute in TE buffer (10 – 20 times) before loading.

We recommend loading 1 mkg of the DNA Ladder per line.



Introduction:

The Hind III digest of lambda DNA yields 8 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | Base Pairs | Base Pairs | Base Pairs |
|------------|------------|------------|------------|
| 23130 | 6557 | 2322 | 564 |
| 9416 | 4361 | 2027 | 125 |

